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## Brochure

This Brochure is additional and is not intended as a substitute for the Guarantee Booklet.

**Manufacturer: B&T S.p.A. via Due Ponti 9, Forlì – Italy**

**Intended use of the product:** Anti-decubitus mattress intended to optimise pressures during rest in order to reduce the risks of prolonged capillary occlusion and decubitus ulcer in individuals with reduced mobility.

Dorelan pressure-reducing mattresses in composite foam (hereinafter referred to simply as mattresses) are suitable to optimise pressures during rest in order to reduce the risks of prolonged capillary occlusion in individuals with reduced mobility. This characteristic is particularly appreciated in the health sector since it plays an important role in preventing pressure ulcers of stage one, if combined with specific changes of position. The reduction of pressure in the support zones is obtained by avoiding a steady posture which excessively stimulates sensitive body zones, if necessary by using accessories (e.g.: pillows).

The product complies with the general safety and performance requirements specified in the Annex I of Regulation (EU) 2017/745.

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### INSTRUCTIONS ON THE CORRECT USE OF THE PRODUCT

It is recommended that the bed base slats are inserted into the frame by means of suitable holders, which allow the slats to be flexible and adjustable, thus providing flexibility to the whole base, including the perimeter.

It is recommended to check the adequacy of the product (i.e. mattress height) by placing a hand between the mattress and the slatted base, in order to check if the body touches the base (unsuitable support) or if it doesn't (suitable support).

Protect the mattress with disposable drawsheets which are air-permeable but impermeable to fluids. The mattress is not fire-resistant. The mattress is not sterile and cannot be sterilized. Do not expose to direct sunlight. Do not spill liquids on the product.

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B&T S.p.A. as the manufacturer and as the only responsible for medical devices placed on the market in its own name, declares, under its own responsibility, that the device to which this declaration refers complies with the general safety and performance requirements specified in the Annex I of the Regulation (EU) 2017/745.

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**Any modification made to the Medical Device unless authorized by B&T Spa cancels the validity of this declaration.**

We hereby declare that this device is sold in a non-sterile package.

We declare that the devices in question comply with the following technical standards:

- EN ISO 14971
- EN ISO 15223-1
- EN 62366

We declare that the undersigned will keep the documentation referred to in Annex I, II and III of Regulation (EU) 2017/745 available to the Competent Authority for a period of 10 years starting from the last production date of the devices in question.

**Forlì - ITALY, 26.5.2021**  
**B&T S.p.A. - Legale Rappresentante**  
**Tura Diano**

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